

Catalog Number: AAA10909

NuPAGE 4-12 % Bis-Tris Gel from Invitrogen

1 mm thick gel


MOPS running buffer

Sample size = 12 ul/lane

See Blue Plus 2 Prestained Standards from Invitrogen

200 Volts constant for 50 min

Stained with Coomassie Brilliant Blue R-250



Lane # 1 2 3 4

Lane #1: 2 μ g C1q, Lot # 622PC1QL + DTT

Lane #2: 4 μ g C1q, Lot # 622PC1QL + DTT

Lane #3: 6 μ g C1q, Lot # 622PC1QL + DTT

Lane #4: See Blue Plus 2 Prestained Standards