

Human Zinc transporter 8 (SLC30A8) autoantibody ELISA kit

PRECISION

Intra-assay Precision (Precision within an assay): CV%<15%

Three samples of known concentration were tested twenty times on one plate to assess.

Inter-assay Precision (Precision between assays):CV%<20%

Three samples of known concentration were tested in twenty assays to assess.

TYPICAL DATA

Test parameter	specification	test result
Positive control	≥0.45	1.195
Negative control	≤0.1	0.037
Positive rate	10, Positive	100%
Negative rate	10, Negative	100%